
GUIDION HYDRO
GUIDION SHORT
FLOWGUIDE

English [EN] Instructions for use - Guide extension family.....	1
Español [ES] Instrucciones de uso - Familia de extensión con guía.....	2
Français [FR] Mode d'emploi - Famille Extension de guide.....	4
Deutsch [DE] Gebrauchsanweisung - Führungserweiterungsreihe.....	5
Italiano [IT] Istruzioni per l'uso - Famiglia Guide Extension.....	6
Türkçe [TR] Kullanım talimatları - K-lavuz uzatma kateteri grubu.....	7
Slovensky [SK] Návod na použitie - Rad predlžovacích vodičiacich katétov.....	9
Português [PT] Instruções de utilização - Família Guide Extension.....	10
Polski [PL] Instrukcja użytkowania - Gama cewników stanowiących przedłużenie prowadnika.....	11
Česky [CS] Návod k použití - řada vodičích nástavců.....	12
Dansk [DA] Brugervejledning - Guideforlængerserien.....	13
Suomi [FI] Käyttöohjeet - Ohjainlaajennustutkitaiden perhe.....	15
Norsk [NO] Bruksanvisning - Ledeforlængersesamilien.....	16
Svenska [SV] Bruksanvisning - Ledarförlångningsfamilj.....	17
[BG] -.....	18
Eesti [ET] Kasutusjuhend - Juhtimispidenduse perekond.....	20
Magyar [HU] Használati utasítás - Vezetőemeghosszabbító család.....	21
Latviešu [LV] Lietošanas pamcība - Vadītāja gaisa pagarinājuma sērija.....	22
Lietuvių [LT] Naudojimo instrukcijos - Kreipiamoji kateterio ilginčių serija.....	23
Română [RO] Instrucțiuni de utilizare - Gama extensibilor de ghidare.....	24
Slovenščina [SL] Navodila za uporabo - Družina s podaljški vodil.....	26

Device group	Device Name	Model	UDI DI
Guidion Rapid Exchange Guide Extension	Guidion 5F	G50F25150	08718481740200
	Guidion 6F	G60F25150	08718481740149
	Guidion 7F	G70F25150	08718481740217
	Guidion 8F	G80F25150	08718481740224
Guidion Short Rapid Exchange Guide Extension	Guidion Short 6F	GS60F	08718481740330
FlowGuide Rapid Exchange Guide Extension	FlowGuide 6F	FG60F	08718481740347

Note: Not all device models are available in every country.

1. [EN] Shaft [ES] Eje [FR] Gaine [DE] Schaft [IT] Asta [TR] Mil [SK] Telo [PT] Haste [PL] Trzon [CS] Důk [DA] Skæft [FI] Varsi [NO] Skæft [SV] Skaf [BG] [ET] Vars [HU] Szár [LV] Vrpsta [LT] Velenas [RO] Tijă [SL] Gred
2. [EN] Push rod [ES] Varilla de empuje [FR] Tige poussoir [DE] Stößel [IT] Pistoncino di spinta [TR] İtme çubuğu [SK] Ťlačidlo [PT] Haste impulsora [PL] Popychacz [CS] Tlačítko [DA] Stadstang [FI] Työntösauna [NO] Skyvestang [SV] Stötstang [BG] [ET] Surumisvarras [HU] Nyomórúd [LV] B distien tīs [LT] St. miklis [RO] Tijă de împingere [SL] Potisni nosilec
3. [EN] Proximal end stop [ES] Final proximal [FR] Butée d'extrémité proximale [DE] Proximaler Endanschlag [IT] Fine corsa prossimale [PT] Proximal uc [SK] Koncovka na proximálné [Basti] [PT] Extremitade proximal [PL] Proksymalna końcówka zatrzymująca [CS] Proximální koncová zarážka [DA] Proksimalt endestop [FI] Proksimaalinen pää [NO] Proksimalt endestopp [SV] Proximalt slutstykke [BG] [ET] Proksimaalne otsak [HU] Proximális csúcs [LV] Proksim. lais gais aizturis [LT] Proksimalin. galin atrama [RO] Capăt terminal proximal [SL] Proksimalni končni del
4. [EN] Soft distal end [ES] Extremo distal blando [FR] Extrémité distale souple [DE] Weiches distales Ende [IT] Estremità distale morbida [PT] Yumulaık distal uc [SK] Mäkký distálny koniec [PT] Extremitade distal maleável [PL] Miękki koniec dystalny [CS] Měkký distální konec [DA] Blød distal ende [FI] Pehmeä distaalinen kärki [NO] Myk distal ende [SV] Mjuk distal änd [BG] [ET] Pehme distaalne ots [HU] Puha distális vég [LV] M ksts dist. lais gais [LT] Minkštas distalinis galas [RO] Capăt distal moale [SL] Mehki distalni konec
5. [EN] Radiopaque a-traumatic tip [ES] Punta distal a-traumatica radiopaca [FR] Extrémité atraumatique radio-opaque [DE] Röntgendichtende atraumatische Spitze [IT] Punta atraumatica radiopaca [TR] Radyopak atravmatik uc [SK] RTG-kontrastný atraumatický hrot [PT] Ponta atraumática radiopaca [PL] Końcówka

radiocieniuja a-traumaticzna [CS] Radiopakní atraumatický hrot [DA] Røntgenfast atraumatisk spids [FI] Radio-opaakki säästävää kärki [NO] Røntgenlett atraumatisk spiss [SV] Röntgentät atraumatisk spet [BG] [ET] Röntgenkontrastne atraumaatiline
otsak [HU] Sugárfogó a-traumás hegy [LV] Røntgenkontrastans, atraumatisk gais [LT] Røntgeno spinduliuose matomas atraumatinis antgalis [RO] Vâr radiopac atraumatic [SL] Radiopaka a-traumatska konica
6. [EN] Exit Marker [ES] Marcador de salida [FR] Marqueur de sortie [DE] Ausgangsmarkierung [IT] Indicatore di uscita [TR] Çıkış Belirteci [SK] Výstupná značka [PT] Marcação de saída [PL] Znacznik wyjścia [CS] Výstupní značka [DA] Udgangsmarker [FI] Poistumismerkki [NO] Utgangsmarker [SV] Utgångsmärken [BG] [ET] Väljumsimärgistus [HU] Kimeneti marker [LV] Izejas marieris [LT] Išėjimo žyma [RO] Indicator de ieşire [SL] Oznajevalec izhoda
7. [EN] For device model FG60F: Side holes [ES] Para FG60F: Agujeros laterales [FR] Pour FG60F: Troux latéraux [DE] Für FG60F: Seitenöffnungen [IT] Per la FG60F: Fori laterali [TR] FG60F için: Yan delikler [SK] Pre katéter FG60F: Bočné otvory [PT] Para o FG60F: Orifícios de deslizeamento [PL] Dotyczy FG60F: Otwory boczne [CS] Pro FG60F: Boční otvory [DA] For FG60F: Sidehuller [FI] FG60F: Sivuaukot [NO] Til FG60F: Sidehull [SV] För FG60F: Sidehål [BG] FG60F: [ET] FG60F-i jaoks Külgmised avad [HU] A FG60F-esetében: Oldalsó furatok [LV] Paredz tas FG60F: Snu atveres [LT] „FG60F” atveju: Šonin s angos [RO] Pentru FG60F: Orificii laterale [SL] Za FG60F: Stranske odprtine
8. [EN] For device models FG60F and GS60F: Radiopaque marker [ES] Para FG60F y GS60F: Marcador radiopaco [FR] Pour FG60F et GS60F: Marqueur radio-opaque [DE] Für FG60F und GS60F: Röntgendichte Markierung [IT] Per la FG60F e il GS60F: Marker radiopaco [TR] FG60F ve GS60F için: Radyopak işaretleyici [SK] Pre katéter FG60F a GS60F: RTG-kontrastná značka [PT] Para o FG60F e o GS60F: Marcador radiopaco [PL] Dotyczy FG60F i GS60F: Znacznik nieprzepuszczający promieni rentgenowskich [CS] Pro FG60F a GS60F: Radiopakní značka [DA] For FG60F og GS60F: Røntgenfast marker [FI] FG60F ja GS60F: Radio-opaakki merkki [NO] Til FG60F og GS60F: Røntgenlett marker [SV] För FG60F och Guidion kort: Röntgentät märkar [BG] FG60F GS60F: [ET] FG60F-i ja GS60F-i jaoks Røntgenkontrastne määrgistus [HU] A FG60F és a GS60F esetében: Sugárfogó marker [LV] Paredz tas FG60F un Guidion Short: Røntgenkontrastans marieris [LT] „FG60F” ir „GS60F” atveju: Røntgeno spinduliuose matoma žyma [RO] Pentru FG60F i GS60F: Marcar radiopac [SL] Za FG60F in GS60F: Radiopaka označevalec

English [EN] Instructions for use - Guide extension family

Device Description

The Guide extension family catheter is a single lumen rapid exchange catheter offered in sizes being compatible with 5F, 6F, 7F and 8F guide catheters and may be placed over either an exchange length or 180cm guide wire. The 150cm long device has a stainless-steel shaft section. The stainless-steel shaft is followed distally by a lumen section. The Guide extension family catheter has a radiopaque distal end which enables visibility while using standard fluoroscopic methods. The device has two positioning marks located at 95cm and 105cm from the distal tip, respectively. The Guide extension family catheter is delivered through a guiding catheter resulting in an inner diameter that is approximately 1 French smaller than the guide catheter. The Guide extension family catheter has a proximal hub which indicates guide catheter compatibility.

STERILE. Sterilized with Gamma Radiation. Non-pyrogenic. For one time use only. Do not re-sterilize.

Eudamed and SSCP

A Summary of the Safety and Clinical Performance (SSCP) for this device is available at <https://ec.europa.eu/tools/eudamed>. Search for the device using the Basic UDI-DI: 871848174GCE0001F8.

The Guide extension family is available in the following size

Model	Compatible guide catheter	Min ID	Guide extension tip OD	Extension shaft
G50F25150	a5F (a0.056" / 1.42mm)	1.04mm	1.34mm	25cm
G60F25150	a6F (a0.070" / 1.78mm)	1.42mm	1.70mm	25cm
G70F25150	a7F (a0.078" / 1.98mm)	1.57mm	1.90mm	25cm
G80F25150	a8F (a0.088" / 2.24mm)	1.80mm	2.16mm	25cm
GS60F, FG60F	a6F (a0.070" / 1.78mm)	1.42mm	1.72mm	15cm

Intended purpose

The Guide extension family catheters are intended to be used in conjunction with guide catheters to access discrete regions of the coronary and/or peripheral vasculature, and to facilitate placement and exchange of guide wires and other interventional devices.

Indications for use

The Guide extension family catheters are indicated for patients suffering from coronary or peripheral artery disease.

Contraindications

The Guide extension family catheters are contraindicated in vessels less than 2.5mm in diameter, vessels in the neurovasculature and the venous system.

Warnings

- For single use only. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn, may result in patient injury, illness or death.
- Handle with care to avoid product damage.
- Use the catheter before the "Use By" date specified on the package.
- Do not use if the package is open or damaged.
- Care should be taken when handling the distal part of the catheter to prevent damage.
- Before inserting the device into the guiding catheter make sure the appropriate exit marker is noted.
- Never advance the Guide extension family catheter into a vessel without a leading guide wire or without confirming location using fluoroscopic guidance. Vessel dissection or perforation may result.
- Never advance the Guide extension family catheter into a vessel with an effective diameter less than 2.5mm. Vessel injury, ischemia, and/or occlusion may result. If pressure in a vessel dampens after inserting the Guide extension family catheter, withdraw the Guide extension family catheter immediately.
- Due to the size and non-tapered tip of the Guide extension family, extreme care must be taken to avoid vessel occlusion and damage to the wall of the vessels through which this catheter passes.
- If resistance is felt during manipulation, determine the cause of the resistance before proceeding.
- Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guide wire against resistance may result in separation of the catheter or guide wire tip, damage to the catheter, or vessel perforation.
- Do not use the Guide extension family during an MRI procedure.

Precautions

- The catheter should be used by physicians trained on the procedures for which the device is intended. The techniques and procedure described do not represent ALL medically accepted protocols, nor are they intended as a substitute for the physician's experience and judgement in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.
- Inspect the Guide extension family catheter prior to use for any bends or kinks. Do not use a damaged catheter. Vessel damage and/or inability to advance or withdraw the catheter may occur.
- The catheter lumen should be thoroughly flushed with heparinized saline prior to use to prevent clot formation.
- Precautions to prevent or reduce clotting should be taken when any catheter is used in the vascular system. Use of systemic heparinization and heparinized sterile solution should be considered.
- Exercise care in handling of the catheter during a procedure to reduce the possibility of accidental breakage, bending or kinking.
- When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response.
- Never advance the Guide extension family catheter more than 15 cm beyond the tip of the guide catheter as the Guide extension family catheter may become lodged in the guide catheter making it difficult to remove.
- The 8F Guidion catheter may be advanced up to 15cm in the peripheral vasculature but should not be advanced more than 10cm in the coronary vasculature
- Do not withdraw an un-deployed stent back into the Guide extension family when the catheter is in the body. Instead, simultaneously pull both the Guide extension family and un-deployed stent back into the guide and remove together.
- Do not torque the catheter more than one (1) full turn.

Reporting of adverse events

The user and / or patient should report any serious incident that has occurred in relation to this Guide extension family to the manufacturer and the competent authority of state in which the user and / or patient is established.

Potential side effects / Adverse events

Possible side effects / adverse events include, but are not limited to the following

- Air embolism
- Intimal disruption
- Arterial dissection
- perforation of the vessel wall
- Vascular occlusion
- Arterial Thrombosis
- Myocardial infarction
- Arterial spasm

Package contains

One Guide extension family catheter

Other items required but not provided

Interventional devices used for the intended purpose of the Guide extension family are described in this IFU. The use and compatibility of the interventional devices

with the Guide extension family is described in the instructions supplied with the respective devices.

Preparation and inspection procedure

Note: Do not expose the catheter to organic solvents (e.g. alcohol).

Note: Store in a cool, dark and dry place

- Prior to use, carefully inspect the Guide extension family catheter packaging and components for damage.
- Using sterile technique, transfer the dispenser coil with the Guide extension family catheter into the sterile field.
- Gently remove the catheter from the packaging hoop.
- Thoroughly flush the lumen of the Guide extension family catheter from the distal tip with heparinized saline solution.
- Wet the outside of the catheter with sterile heparinized saline to activate the hydrophilic coating.

Deployment procedure

Deploy the Guide extension family catheter according to the following steps

- Secure the previously inserted guide wire and backload the distal end of the Guide extension family catheter onto the guide wire and advance until it is just proximal to the hemostasis valve.
- Open the hemostasis valve and advance the Guide extension family catheter through the hemostasis valve and into the guide catheter.
- Gently advance the Guide extension family catheter until the distal end of the Guide extension family catheter is approximately aligned with the tip of the guiding catheter.
- Under fluoroscopy, advance the Guide extension family catheter beyond the distal tip of the guide catheter and into the desired location within the vessel.

Note: When using models with side holes, partial contrast injection through the guiding catheter will flow out through the side holes.

When using models with side holes, blood entering in the coronary where the guide extension is inserted, will partially flow through the side holes to the vessels distally of the guide extension tip

Warning: Never advance the Guide extension family catheter into a vessel with an effective diameter less than 2.5mm. Vessel injury, ischemia, and/or occlusion may result. If pressure in a vessel dampens after inserting the Guide extension family catheter, withdraw the Guide extension family catheter until the pressure returns to normal

Warning: Due to the size and non-tapered distal end of the Guide extension family, extreme care must be taken to avoid vessel occlusion and damage to the wall of the vessels through which this catheter passes.

Note: The exit markers on the body of the shaft of the catheter can be used to estimate the insertion depth.

- Using fluoroscopy, confirm the desired position of the Guide extension family catheter in the vessel.
 - If performing an interventional procedure, backload the interventional device over the in place guide wire and advance the device through the guide catheter and Guide extension family catheter into the desired vascular space.
- Note:** If a second wire is used during the intervention and encounters resistance within the guide catheter, pull the wire back several centimeters and slowly re-advance.
- Tighten the Y-adaptor hemostasis valve securely on the proximal shaft of the Guide extension family catheter to prevent back-bleeding.
 - Perform the catheterization procedure. After completing the procedure, remove the Guide extension family catheter prior to removing the guide catheter from the vessel.

Guide extension family catheter removal

- Keep guide wire in position
- Gently retract the catheter out of the guiding catheter and of the guide wire using the standard technique.

Product disposal.

After use, this product may be a potential biohazard. Handle and dispose of it in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

References

The physician should consult recent literature on current medical practice.

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[EN] Disclaimer of warranty and limitation of Remedy

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the IMDS product(s) described in this publication. Under no circumstances shall IMDS be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind IMDS to any representation or warranty except as specifically set forth herein.

Description or specifications in IMDS printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. IMDS will not be responsible for any direct, incidental or consequential damages resulting from the reuse of the product.

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