RECROSS

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English [EN] Instructions for use - ReCross OTW Dual Lumen Guide Wire Support Catheter

Device Description

The ReCross is 0.014" Over the Wire (OTW) dual lumen guide wire support catheter with two lumens that run the length of the catheter. The ReCross is intended to support guide wires and allows for exchange of guide wires. The ReCross provides a conduit for local delivery of diagnostic agents. The ReCross has an effective length of 140 cm and is compatible with 0.014"/0.36mm or smaller guide wires. The ReCross has two depth markings located at 95 cm and 105 cm from the distal tip. The ReCross has a radiopaque marker identifying the distal end of the longest lumen (tip lumen) of the catheter, a second radiopaque marker located 8mm from the distal tip identifies the distal end of the shorter lumen (stylet lumen) of the catheter, and a third radiopaque marker located 12mm from the distal end locating an additional exit port in the tip lumen 180° opposite of the end of the stylet lumen. The distal 25 cm of the shaft is fitted with a hydrophilic coating. The entire lumens of the ReCross are coated with MDX (silicone) coating. To ease catheter loading into the hemostasis device and the guiding catheter a removable stylet is placed in the Stylet lumen. The ReCross is provided with two extensions at the proximal end with indication of the lumen on the luer connections.

STERILE. Sterilized with Gamma Radiation. Non-pyrogenic. For one time use only. Do not re-sterilize.

Eudamed and SSCP

A Summary of the Safety and Clinical Performance (SSCP) for this device is available at <https://ec.europa.eu/tools/eudamed>. Search for the device using the Basic UDI-DI: 871848174GWS0003UA.

The ReCross is available in the following size

Model	Compatible guide catheter	Max Guide wire size	ReCross ID
RC1414025	æbF (æd.056" /1.42mm)	0.36mm (0.014")	0.40mm (0.016")

Intended purpose

The ReCross is a dual lumen OTW guide wire support catheter intended to support a guide wire during access of coronary and/or peripheral vasculature, allows for exchange of guide wires and provides a conduit for the delivery of diagnostic contrast agents.

Indications for use

The ReCross is indicated for patients suffering from coronary or peripheral artery disease.

Contraindications

The catheter is not indicated for use in the neurovasculature.

Warnings

- For single use only. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn, may result in patient injury, illness or death.
- Use the catheter before the "Use By" date specified on the package.
- Handle with care to avoid product damage.
- Do not use the ReCross catheter if the inner package is open or damaged.
- This device is coated with a hydrophilic coating at the distal end of the device for a length of 25 cm and is coated with a silicone coating throughout both lumens of the catheter. Please refer to Preparation and inspection procedure in this IFU for further instruction how to prepare and use this device to ensure it performs as intended. Failure to abide by the warnings in this IFU might result in damage to the device coating, which may necessitate intervention or results in serious adverse events / side effects.
- Do not advance the ReCross catheter without a guide wire in place through the distal tip of the tip lumen. Advancement of the catheter without a guide wire in the distal tip may result in vessel damage.
- Before inserting the device into the guiding catheter make sure the appropriate exit marker is noted.
- Insert / advance the device into the vasculature system with care.
- To reduce the potential for vessel damage the support catheter should not be inserted in vessels with an inner diameter smaller than the outer diameter of the support catheter.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation.
- Do not use a power injection system when injecting through the ReCross because this might cause catheter damage or vessel damage.

- Care should be taken when handling the distal part of the catheter to prevent damage.
- At guide wire removal from OTW lumen ensure catheter is secured in position.
- Do not advance the ReCross catheter into cerebral vasculature. Advancement in the cerebral vasculature might lead to stroke and/or cerebral infarct, which may result in patient injury, illness or death.
- If an appropriately-sized diameter guidewire cannot be passed through the catheter, do not attempt to resolve the blockage by flushing the catheter in vivo. Catheter rupture and arterial injury may result. Identify and resolve the cause of the blockage or replace the catheter with a new one.
- Never advance, withdraw, or rotate an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter against resistance may result in separation of the catheter or guidewire tip, another device damage, or vessel injury
- The ReCross it is not intended to be rotated more than one turn (360°) to prevent distal guide wire entanglement.

Precautions

- Prior to use, the catheter should be examined to verify functionality and ensure that its size is suitable for the specific procedure for which it is to be used.
- Appropriate anticoagulant / antiplatelet therapy should be used during this procedure.
- During a procedure the catheter should be handled with care to reduce the possibility of accidental breakage, bending or kinking.
- Care should be taken not to crush the catheter. Excessive tightening of the hemostatic valve onto the catheter shaft may result in damage to the guidewire lumen and difficulty when inserting the catheter or guidewires.
- The catheter should be used by physicians trained on the procedures for which the device is intended. The techniques and procedure described do not represent ALL medically accepted protocols, nor are they intended as a substitute for the physician's experience and judgement in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

Reporting of adverse events

The user and / or patient should report any serious incident that has occurred in relation to this ReCross to the manufacturer and the competent authority of state in which the user and / or patient is established.

Potential side effects / Adverse events

Possible side effects / adverse events include, but are not limited to the following

- Death
- Acute myocardial infarction
- Total occlusion of the coronary artery or bypass graft.
- Coronary vessel dissection, perforation, rupture or injury
- Restenosis of the dilated vessel.
- Hemorrhage or hematoma.
- Unstable angina.
- Arrhythmias, including ventricular fibrillation.
- Drug reactions, allergic reaction to contrast medium.
- Hypo / hypertension.
- Infection
- Coronary artery spasm
- Arteriovenous fistula.
- Pulmonary embolism
- Pulmonary infarct
- Myocardial embolism
- Embolitic stroke
- Cerebral infarct
- Tissue necrosis
- Inflammation
- Thrombosis
- Local or systemic injection
- Catheter fraction with tip separation and distal embolization

Contents

One (1) ReCross OTW Guide Wire Support Catheter. The Stylet is used as an accessory to ease catheter loading into the hemostasis device and guide catheter

Other items required but not provided

Interventional devices used for the intended purpose of the ReCross are described in this IFU. The use and compatibility of the interventional devices with the ReCross is described in the instructions supplied with the respective devices.

Preparation and inspection procedure

Note: Do not expose the catheter to organic solvents (e.g. alcohol).

Note: Store in a cool, dark and dry place

Warning: Do not use if the inner package is open or damaged.

- Gently remove the catheter stored in hoop from the carton and pouch.
- Gently detach the stylet from the hub and remove the stylet from the stylet lumen.
- Gently detach the hub of the catheter from the packaging clip.
- Attach a syringe filled with a minimum of 1 cc of sterile heparinized saline to the stylet hub.
- Flush the stylet lumen of the catheter.

- Remove the syringe
- Attach a syringe filled with a minimum of 1 cc of sterile heparinized saline to the tip hub.
- Flush the tip lumen of the catheter.
- Remove the syringe
- Gently remove the catheter from the packaging hoop.
- Inspect the catheter on kinks or other damages.
- Wet the outside of the catheter with sterile heparinized saline to activate the hydrophilic coating.

Assembly and Insertion Procedure

- Make sure the guide wire is soaked in heparinized saline.
- Backload the guide wire that is already in place in the vasculature into the tip lumen of the catheter.

Note: To support in loading the catheter the stylet or a second guide wire can be inserted in the stylet lumen.

- Open the hemostasis device mounted on the guiding catheter

Note: Use a 5F Guiding Catheter as a minimum.

Warning: Make sure the appropriate exit marker is noted.

- Gently advance the catheter over the guide wire into the guiding catheter until the tips of the guiding catheter and the catheter are approximately aligned.

Warning: Do not advance the ReCross catheter without a guide wire in place through the distal tip of the tip lumen. Advancement of the catheter without a guide wire in the distal tip may result in vessel damage.

Note: The exit markers on the body of the shaft of the catheter can be used to estimate the insertion depth.

- Carefully advance the catheter under fluoroscopy until the lesion.

Warning: If resistance is felt during manipulation, determine the cause of the resistance before proceeding.

- Carefully position the distal tip or one of the exit ports at the required position

Guide wire assembly and (re-)introduction in a lumen.

- Make sure the guide wire is soaked in heparinized saline.
- Gently introduce the guide wire into the tip hub or stylet hub of the catheter.
- Carefully advance the guide wire up to the distal part of the catheter and the lesion.

Warning: If resistance is felt during manipulation, determine the cause of the resistance before proceeding.

- Carefully manipulate the guide wire in the selected arteries.

Warning: If resistance is felt during manipulation, determine the cause of the resistance before proceeding.

Guide wire removal from a lumen

- Keep the catheter in its position.
- Remove the guide wire carefully.

Administration of contrast medium

- Position the ReCross with the distal tip or one of the exit ports at the place of interest.
- Keep ReCross in its position.
- Remove the guide wire from the catheter.
- Attach an appropriate syringe with contrast medium to tip lumen hub or stylet lumen hub.
- Inject the contrast medium.
- Remove the syringe
- Re-insert the guide wire.

Catheter removal

- Secure both guide wires using standard guide wire exchange techniques.
- Keep the guide wires in position.
- Open the hemostasis device mounted on the guiding catheter
- Carefully retract the catheter out of the guiding catheter and of the guide wire using the standard technique.

Product disposal.

After use, this product may be a potential biohazard. Handle and dispose of it in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

References

The physician should consult recent literature on current medical practice.

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[EN] Disclaimer of warranty and limitation of Remedy

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the IMDS product(s) described in this publication. Under no circumstances shall IMDS be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind IMDS to any representation or warranty except as specifically set forth herein.

Description or specifications in IMDS printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. IMDS will not be responsible for any direct, incidental or consequential damages resulting from the reuse of the product.

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IMDS Operations B.V.