



CerebrAXTM

Stent Retriever



HEMO BIOENGINEERING PTE. LTD.

Website: www.hemocorp.com E-mail: enquire@hemocorp

2024 HEMO SINGAPORE. All rights reserved. 2024-G2-V1.00

Longer effective Length

More favourable clot-stent engagement thereby increasing probability of successful clot retrieval


Comprehensive OD size range

Greater size selection facilitates selection of different stent sizes and lengths for varying vessel diameters

Midpoint radiopaque markers

Physicians will be able to judge the relative position of the stent with respect to the thrombus, thereby facilitating the placement of stent.

Product Structure



- | | |
|-----------------------------------|--|
| 1. Distal radiopaque markers | A. Nominal diameter of stent |
| 2. Middle radiopaque markers | B. Distance between middle and distal radiopaque markers |
| 3. Self-expanding stent (nitinol) | C. Effective length of self-expanding stent |
| 4. Proximal radiopaque markers | D. Distance between proximal and distal radiopaque markers |
| 5. Introducer sheath (PTFE) | E. Total length of stent retriever |
| 6. Push wire (nitinol) | |

Product Specification

Model and Specification	Dimension of Stent (mm)	Nominal Diameter of Stent (mm)	Compatible Vessel OD (mm)	Effective Length of Stent (mm)	Total Length (mm)	No. of Radiopaque Markers			Min ID of Microcatheter (mm/inch)
						Proximal	Middle	Distal	
SARD-3-10	3×10	3	2.0-3.0	10	1850	1	2	2	0.53/0.021
SARD-3-20	3×20	3	2.0-3.0	20	1850	1	2	2	0.53/0.021
SARD-4-20	4×20	4	2.5-4.0	20	1850	1	3	3	0.53/0.021
SARD-4-40	4×40	4	2.5-4.0	40	1850	1	3	3	0.53/0.021
SARD-5-20	5×20	5	3.5-4.5	20	1850	1	3	4	0.69/0.027
SARD-5-30	5×30	5	3.5-4.5	30	1850	1	3	4	0.69/0.027
SARD-6-20	6×20	6	4.0-5.5	20	1850	1	3	4	0.69/0.027
SARD-6-30	6×30	6	4.0-5.5	30	1850	1	3	4	0.69/0.027